


Waves Analyser


FlameinSounds Waves Analyser
was originally designed to display a relationship between
Heart Rate Variations
&
the **beat** generated in the Binaural sound wave.

When you are very relax and centered - a state
inducing a much slow and much deep kind of breathing,
you body can enter in resonance with the beat (rhythm)
generated by the Binaural sound waves.


The Waves Analyser can also be used
without the sound generator.
In this case, only heart waves will be displayed.




Waves Analyser




Display Options




Select the curves you want to display/hide.



Allows to select a predefined Time Scale of the display. This Time Scale is saved for the next use of the app.




Swipe RIGHT/LEFT
the middle of the screen
to change the Time Scale
during a live session or a playback.




1 min Time Scale

With HRV sensors,
it is possible to display
Heart Rate Variations
&
beat generated
in the Binaural sound
wave.



5 min Time Scale
of the same recording



With PPG sensor(BERRY),
it is possible to display
Heart Rate Variations
&
beat generated
in the Binaural sound
wave.
&
Blood Pressure Pulses
&
Mayer Wave
(Sacro-Cranial Wave)

